



**Name:** [Anonymous]

**Product:** NANI™, Telepathy™ NBS, Telepathy™ CCHD, Telepathy™ EHDI, OZ eScreener Plus™ (eSP™)

This testimonial is from a Texas Director of Women's and Children's Services, and like many of our hospital customers, the testimonial could only be provided anonymously. She is passionate about enhancing workflow and improving patient safety; therefore, she recommends the full OZ solution.



OZ solutions reduce human errors, improve patient safety and reduce the time needed to perform screenings. Prior to OZ NANI, we looked in multiple places to obtain the information needed to complete the specimen card. NANI pulls the information into a form, saving time and increasing accuracy. Telepathy NBS barcoding of the specimen ID has eliminated the errors that occur when information is entered manually. Telepathy CCHD has improved protocol accuracy because it standardizes the process and requires obtaining multiple values. It is very easy to perform and much more accurate. Recently a newborn did not pass the CCHD screening. A follow-up cardiology appointment was scheduled prior to discharge. We had piece of mind knowing we provided the patient and their family the best care. OZ Solutions absolutely improves workflow and patient care.



In approximately thirty minutes we perform four screenings: hearing, CCHD, blood spot and transcutaneous bilirubin measurement. We have technicians who perform the screenings and then they report to the nurses. OZ NANI™ initiates the patient record and is utilized for all three screenings. Our nursing staff loves being able to see a complete view of each baby's newborn screening results in one location: hearing, CCHD and blood spot screening.



**oz systems**

Screening Solutions.  
Improving New Lives.

We are committed to delivering the best standards-based technologies to save more lives and help make the healthcare experience caring, smart and trust-worthy.

888.727.3366 | [www.ozsystems.com](http://www.ozsystems.com)



@OZSystems



/OZSystems



/OZSystems